


Limb holder

Applications & Indications


Operating rooms, hospital wards, ICU, emergency rooms. The limb holder is designed to limit and control the movements of confused patients.



The Velcro closure allows an easy application and quick release of the device. The device must be applied on the upper (wrists) or lower (ankles) limbs only. Washable.

2100	Fast opening limb holder – adult – standard	60 pcs/box
2100TM	Fast opening limb holder – adult breathable, in medical TNT	60 pcs/box
2150	Fast opening limb holder – paediatric	60 pcs/box

2200	Reinforced limb holder with ring – standard	60 pcs/box
2200TM	Reinforced limb holder with ring breathable, in medical TNT	60 pcs/box



Limb holder - glove

Applications & Indications

The device is designed with individual finger slots to limit finger and hand movement. Washable.

2300	Limb holder with protection glove	16 pcs/box
------	-----------------------------------	------------

Elbow – heel protector

Applications & Indications

The device is designed to prevent pressure sores and must be positioned on the elbows and on the heels. Washable.

4100	Elbow-heel protector – standard	36 pcs/box
4100S	Elbow-heel protector – extra padding breathable, in medical TNT	36 pcs/box



4200	Open elbow-heel protector – long	20 pcs/box
------	----------------------------------	------------

LATEX FREE

Class I Medical Device

Made in Italy

Copyright 2003 – COREMEC s.r.l. Edition SP40/r1/eng/2006 - Italy

COREMEC s.r.l.

Manufacturing Plant & Administration
Via Cacciatori, 74/76 - 10042 Nichelino (TO) - Italy
Phone ++39 011.35.89.481
Fax ++39 011.35.89.483
www.coremec.it e-mail coremec@coremec.it



UNI EN ISO 9001:2000 Certified

COREMEC s.r.l.
Export Department

e-mail: export@coremec.it