

Intensive Care Unit

Neonatology

Home Care

Fixing of Nasal Probes

The use of nasal and gastric probes is a frequent procedure in hospitals. The probe runs into the patient's stomach through their nose. Therefore, the probe has to be safely secured onto the patient's nasal bridge. Our Transafix Nasal Fixing System offers a perfect solution.

The set consists of two parts, the nasal clip, available in skin colored fleece or hydrocolloid, and the nasal strip.



The components are connected via a secure Velcro fastener. This way, when exchanging the nasal or gastric probe, only the strip has to be removed; the nasal clip can stay on the nose for up to 3-5 days. This reduces skin irritations and infections and prevents a complex follow-up care.

While both materials, the fleece and hydrocolloid, are very well tolerated by the skin, hydrocolloid is especially recommended for sensitive skin.



The clips and strips are available in a handy dispenser of 20 pieces each. Therefore, the product can be adjusted individually to the needs of each patient.

Benefits

- Easy handling
- Adjustable to the patient's needs
- Handy dispenser
- Easy removal with no residue
- Maximum patient comfort
- Gentle on sensitive skin
- High quality

Nasal Fixing System

Transafix, 1 Nasal Clip - hydrocolloid, non-sterile

Article-No. 8210100

Transafix, 1 Nasal Clip - fleece, non-sterile

Article-No. 8210200

Transafix, 1 Nasal Strip, non-sterile

Article-No. 8210300

Fixing of Nasal Probes

Overview

Intensive Care Unit

Neonatology

Home Care

- Easy connecting and disconnecting of the nasal clip and strip via Velcro fastener
- The nasal clip is not removed when exchanging the probe, only the nasal strip
- The probe is free in the nose, preventing decubitus on the nasal septum

Field of application:

Gastric feeding tubes, oxygen probes, etc.

Application duration:Fleece approx. 3 days, hydrocolloid approx. 5 days,
Nasal Strip as required

Securing the nasal strip on the probe






Illustration of the products in original size

**Transafix, 1 Nasal Clip - hydrocolloid, sterile**

Article-No. 8210100

Transafix, 1 Nasal Clip - fleece, non-sterile

Article-No. 8210200

**Transafix, 1 Nasal Strip, non-sterile**

Article-No. 8210300