

## Belt for lifting and transfer patients

**Ref. CST - M**

**Ref. CST - L**



This belt for lifting and transferring patients represents an effective, comfortable and safe solution for positioning the patient, with limited autonomy, from the patient's bed to the wheelchair and vice-versa. With its 2 side handles and 2 posterior handles, it allows the operators to move the safely the patient, avoiding the risk of sudden movements. Thanks to its particular configuration, the belt can be worn from the operator and from the patient. It is recommended as support also during walking exercises. The internal side of the belt is manufactured in non-slip tissue. The closure is guaranteed through a plastic trident buckle.



### Features:

#### Dimensions:

- CST M : cm 13 x 102 adjustable from cm 74 to 114
- CST L : cm 13 x 102 adjustable from cm 89 to 140

#### Washing instructions :

- Washable at max 70°C with water and mild soap
- Dry clean at max 70°C

Medical Device Class I

Copyright 2009 – Edition COREMEC S.r.l. SP64\_ r0\_2010\_All rights reserved